

## Product Datasheet

### Mouse Interferon gamma ELISA Kit (orb407745)

|                          |                                                       |
|--------------------------|-------------------------------------------------------|
| <b>Description</b>       | Mouse Allograft inflammatory factor 1(AIF1) ELISA kit |
| <b>Reactivity</b>        | Mouse                                                 |
| <b>Range</b>             | 1.10 ng/ml-800 ng/ml.                                 |
| <b>Target</b>            | AIF1                                                  |
| <b>Note</b>              | For research use only                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                          |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                     |
| <b>Uniprot ID</b>        | <b>070200</b>                                         |
| <b>Sensitivity</b>       | 1.10 ng/ml.                                           |
| <b>Expiration Date</b>   | Please enquire.                                       |