

Product Datasheet

Mouse ANGPTL6 ELISA Kit (orb407759)

Description	Mouse Angiopoietin-related protein 6(ANGPTL6) ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	ANGPTL6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8R0Z6
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.