

Product Datasheet

Mouse ANGPTL7 ELISA Kit (orb407760)

Description	Mouse Angiopoietin-related protein 7 (ANGPTL7) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	ANGPTL7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8R1Q3
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.