

Product Datasheet

Mouse PYCARD ELISA Kit (orb407778)

Description	Mouse Apoptosis-associated speck-like protein containing a CARD(PYCARD) ELISA kit
Reactivity	Mouse
Range	7.81 pg/mL-500 pg/mL
Target	PYCARD
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9EPB4
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.