

Product Datasheet

Mouse BTLA ELISA Kit (orb407787)

Description	Mouse B- and T-lymphocyte attenuator(BTLA) ELISA kit
Reactivity	Mouse
Range	0.625 ng/mL-40 ng/mL
Target	BTLA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q7TSA3
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.