

Product Datasheet

Mouse PDGFR beta ELISA Kit (orb407804)

Description	Mouse Beta-type platelet-derived growth factor receptor(PDGFRB) ELISA kit
Reactivity	Mouse
Range	23.44 pg/mL-1500 pg/mL
Target	PDGFRB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P05622
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.