

Product Datasheet

Mouse BMPER ELISA Kit (orb407806)

Description	Mouse BMP-binding endothelial regulator protein(BMPER) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	BMPER
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8CJ69
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.