

Product Datasheet

Mouse PLA2G5 ELISA Kit (orb407830)

Description	Mouse Calcium-dependent phospholipase A2(PLA2G5) ELISA kit
Reactivity	Mouse
Range	62.5 pg/ml-4000 pg/ml
Target	PLA2G5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P97391
Sensitivity	15.6 pg/ml
Expiration Date	Please enquire.