

Product Datasheet

Mouse CLCF1 ELISA Kit (orb407844)

Description	Mouse Cardiotrophin-like cytokine factor 1 (CLCF1) ELISA kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Target	CLCF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9QZM3
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.