

Product Datasheet

Mouse CPT1A ELISA Kit (orb407845)

Description	Mouse Carnitine O-palmitoyltransferase 1, liver isoform(CPT1A) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	CPT1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell lysates, tissue homogenates
Uniprot ID	P97742
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.