

Product Datasheet

Mouse CHRDL1 ELISA Kit (orb407893)

Description	Mouse Chordin-like protein 1(CHRDL1) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	CHRDL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q920C1
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.