

## Product Datasheet

### Mouse Claudin 5 ELISA Kit (orb407899)

|                          |                                                    |
|--------------------------|----------------------------------------------------|
| <b>Description</b>       | Mouse Claudin-5(CLDN5) ELISA kit                   |
| <b>Reactivity</b>        | Mouse                                              |
| <b>Range</b>             | 18.75 pg/mL-1200 pg/mL                             |
| <b>Target</b>            | CLDN5                                              |
| <b>Note</b>              | For research use only                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                       |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates, cell lysates |
| <b>Uniprot ID</b>        | <b>O54942</b>                                      |
| <b>Sensitivity</b>       | 4.68 pg/mL                                         |
| <b>Expiration Date</b>   | Please enquire.                                    |