

Product Datasheet

Mouse C1QTNF3 ELISA Kit (orb407937)

Description	Mouse Complement C1q tumor necrosis factor-related protein 3(C1QTNF3) ELISA kit
Reactivity	Mouse
Range	3.12 ng/mL-200 ng/mL
Target	C1QTNF3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9ES30
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.