

Product Datasheet

Mouse CRHBP ELISA Kit (orb407957)

Description	Mouse Corticotropin-releasing factor-binding protein(CRHBP) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	CRHBP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q60571
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.