

Product Datasheet

Mouse CLEC18A ELISA Kit (orb407961)

Description	Mouse C-type lectin domain family 18 member A(CLEC18A) ELISA kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	CLEC18A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q7TSQ1
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.