

Product Datasheet

Mouse CXCL3 ELISA Kit (orb407967)

Description	Mouse C-X-C motif chemokine 3(CXCL3) ELISA kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	CXCL3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6W5C0
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.