

Product Datasheet

Mouse CYP2E1 ELISA Kit (orb407978)

Description	Mouse Cytochrome P450 2E1(CYP2E1) ELISA kit
Reactivity	Mouse
Range	39 pg/mL-2500 pg/mL
Target	CYP2E1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q05421
Sensitivity	9.75 pg/mL
Expiration Date	Please enquire.