

Product Datasheet

Mouse ADAM17 ELISA Kit (orb408004)

Description	Mouse Disintegrin and metalloproteinase domain-containing protein 17(ADAM17) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	ADAM17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z0F8
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.