

Product Datasheet

Mouse SMURF1 ELISA Kit (orb408024)

Description	Mouse E3 ubiquitin-protein ligase SMURF1(SMURF1) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	SMURF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9CUN6
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.