

Product Datasheet

Mouse ENPP1 ELISA Kit (orb408028)

Description	Mouse Ectonucleotide pyrophosphatase/phosphodiesterase family member 1(ENPP1) ELISA kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	ENPP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P06802
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.