

Product Datasheet

Mouse EFEMP1 ELISA Kit (orb408031)

Description	Mouse EGF-containing fibulin-like extracellular matrix protein 1(EFEMP1) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	EFEMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8BPB5
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.