

Product Datasheet

Mouse EDIL3 ELISA Kit (orb408033)

Description	Mouse EGF-like repeat and discoidin I-like domain-containing protein 3(EDIL3) ELISA kit
Reactivity	Mouse
Range	28 pg/mL-1800 pg/mL
Target	EDIL3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O35474
Sensitivity	7 pg/mL
Expiration Date	Please enquire.