

Product Datasheet

Mouse EGFL6 ELISA Kit (orb408056)

Description	Mouse Epidermal growth factor-like protein 6(EGFL6) ELISA kit
Reactivity	Mouse
Range	93.75 pg/mL-6000 pg/mL
Target	EGFL6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9JJZ5
Sensitivity	23.43 pg/mL
Expiration Date	Please enquire.