

Product Datasheet

Mouse EGFL7 ELISA Kit (orb408057)

Description	Mouse Epidermal growth factor-like protein 7(EGFL7) ELISA kit
Reactivity	Mouse
Range	7.8 ng/mL-500 ng/mL
Target	EGFL7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	Q9QXT5
Sensitivity	1.95 ng/mL
Expiration Date	Please enquire.