

Product Datasheet

Mouse FBXO32 ELISA Kit (orb408072)

Description	Mouse F-box only protein 32(FBXO32) ELISA kit
Reactivity	Mouse
Range	37.5 pg/mL-2400 pg/mL
Target	FBXO32
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9CPU7
Sensitivity	9.375 pg/mL
Expiration Date	Please enquire.