

Product Datasheet

Mouse FGF23 ELISA Kit (orb408082)

Description	Mouse Fibroblast growth factor 23(FGF23) ELISA kit
Reactivity	Mouse
Range	23.44 pg/mL-1500 pg/mL
Target	FGF23
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9EPC2
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.