

Product Datasheet

Mouse FNDC5 ELISA Kit (orb408089)

Description	Mouse Fibronectin type III domain-containing protein 5(FNDC5) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	FNDC5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8K4Z2
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.