

Product Datasheet

Mouse LGALS3BP ELISA Kit (orb408118)

Description	Mouse Galectin-3-binding protein(LGALS3BP) ELISA kit
Reactivity	Mouse
Range	0.39 ng/mL-25 ng/mL
Target	LGALS3BP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q07797
Sensitivity	0.09 ng/mL
Expiration Date	Please enquire.