

Product Datasheet

Mouse GFRAL ELISA Kit (orb408133)

Description	Mouse GDNF family receptor alpha-like(GFRAL) ELISA kit
Reactivity	Mouse
Range	25 pg/mL-1600 pg/mL
Target	GFRAL
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q6SJE0
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.