

Product Datasheet

Mouse SERPINE2 ELISA Kit (orb408135)

Description	Mouse Glia-derived nexin(SERPINE2) ELISA kit Only one 24 Test trial size kit can be applied for each...
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Target	SERPINE2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q07235
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.