

Product Datasheet

Mouse GLP1R ELISA Kit (orb408136)

Description	Mouse Glucagon-like peptide 1 receptor(GLP1R) ELISA kit
Reactivity	Mouse
Range	0.45 ng/ml - 30 ng/ml
Target	GLP1R
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O35659
Sensitivity	0.11 ng/ml
Expiration Date	Please enquire.