

Product Datasheet

Mouse GAD1 ELISA Kit (orb408140)

Description	Mouse Glutamate decarboxylase 1(GAD1) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	GAD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P48318
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.