

Product Datasheet

Mouse PLA2G12A ELISA Kit (orb408164)

Description	Mouse Group XIIA secretory phospholipase A2(PLA2G12A) ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	PLA2G12A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9EPR2
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.