

Product Datasheet

Mouse HSP90AB1 ELISA Kit (orb408180)

Description	Mouse Heat shock protein HSP 90-beta(HSP90AB1) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	HSP90AB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P11499
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.