

Product Datasheet

Mouse HAO1 ELISA Kit (orb408204)

Description	Mouse Hydroxyacid oxidase 1(HAO1) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	HAO1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9WU19
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.