

Product Datasheet

Mouse IGFBP5 ELISA Kit (orb408211)

Description	Mouse Insulin-like growth factor-binding protein 5(IGFBP5) ELISA kit
Reactivity	Mouse
Range	78 pg/ml -5000pg/ml
Target	IGFBP5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q07079
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.