

Product Datasheet

Mouse IGFBP7 ELISA Kit (orb408212)

Description	Mouse Insulin-like growth factor-binding protein 7(IGFBP7) ELISA kit
Reactivity	Mouse
Range	12.5 pg/mL-800 pg/mL
Target	IGFBP7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	Q61581
Sensitivity	3.12 pg/mL
Expiration Date	Please enquire.