

Product Datasheet

Mouse IGFALS ELISA Kit (orb408213)

Description	Mouse Insulin-like growth factor-binding protein complex acid labile subunit(IGFALS) ELISA kit
Reactivity	Mouse
Range	3.12 ng/mL-200 ng/mL
Target	IGFALS
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P70389
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.