

Product Datasheet

Mouse INSL5 ELISA Kit (orb408215)

Description	Mouse Insulin-like peptide INSL5(INSL5) ELISA kit
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Target	INSL5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9WUG6
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.