

Product Datasheet

Mouse ST2 ELISA Kit (orb408240)

Description	Mouse Interleukin-1 receptor-like 1(IL1RL1) ELISA kit
Reactivity	Mouse
Range	78 pg/ml-5000 pg/ml
Target	IL1RL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P14719
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.