

Product Datasheet

Mouse IL10RA ELISA Kit (orb408242)

Description	Mouse Interleukin-10 receptor subunit alpha(IL10RA) ELISA kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	IL10RA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q61727
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.