

## Product Datasheet

### Mouse IL20 ELISA Kit (orb408258)

<b>Description</b>	Mouse Interleukin-20(IL20) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25 pg/ml-2000pg/ml
<b>Target</b>	IL20
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9JKV9</a>
<b>Sensitivity</b>	7.8 pg/ml
<b>Expiration Date</b>	Please enquire.