

Product Datasheet

Mouse IL22RA2 ELISA Kit (orb408261)

Description	Mouse Interleukin-22 receptor subunit alpha-2(IL22RA2) ELISA kit
Reactivity	Mouse
Range	46.88 pg/mL-3000 pg/mL
Target	IL22RA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q80XF5
Sensitivity	11.7 pg/mL
Expiration Date	Please enquire.