


Product Datasheet

Cystatin C antibody (orb40828)

Description	Cystatin C antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Troponin I, fast skeletal muscle protein (1-100AA)
Target	TNNI2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48788
Dilution Range	IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human heart tissue.