

Product Datasheet

Mouse LGI3 ELISA Kit (orb408303)

Description	Mouse Leucine-rich repeat LGI family member 3(LGI3) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	LGI3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8K406
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.