

Product Datasheet

Mouse LAG3 ELISA Kit (orb408327)

Description	Mouse Lymphocyte activation gene 3 protein(LAG3) ELISA kit
Reactivity	Mouse
Range	1.25 ng/mL-80 ng/mL
Target	LAG3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q61790
Sensitivity	0.31 ng/mL
Expiration Date	Please enquire.