

Product Datasheet

Mouse PRCP ELISA Kit (orb408334)

Description	Mouse Lysosomal Pro-X carboxypeptidase(PRCP) ELISA kit
Reactivity	Mouse
Range	37.5 pg/mL-2400 pg/mL
Target	PRCP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q7TMR0
Sensitivity	9.37 pg/mL
Expiration Date	Please enquire.