

Product Datasheet

Mouse MMP12 ELISA Kit (orb408339)

Description	Mouse Macrophage metalloelastase(MMP12) ELISA kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	MMP12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture homogenates
Uniprot ID	P34960
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.