

## Product Datasheet

### Mouse MSR1 ELISA Kit (orb408340)

|                          |                                                                    |
|--------------------------|--------------------------------------------------------------------|
| <b>Description</b>       | Mouse Macrophage scavenger receptor types I and II(MSR1) ELISA kit |
| <b>Reactivity</b>        | Mouse                                                              |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL                                               |
| <b>Target</b>            | MSR1                                                               |
| <b>Note</b>              | For research use only                                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                       |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                  |
| <b>Uniprot ID</b>        | <b>P30204</b>                                                      |
| <b>Sensitivity</b>       | 0.039 ng/mL                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                    |