

Product Datasheet

Mouse UCP1 ELISA Kit (orb408375)

Description	Mouse Mitochondrial brown fat uncoupling protein 1(UCP1) ELISA kit
Reactivity	Mouse
Range	4.68 pg/mL-300 pg/mL
Target	UCP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P12242
Sensitivity	1.17 pg/mL
Expiration Date	Please enquire.